

Monitor Yoke

Medical Display Mount with Swivel

- Standard swivel mount feature that allows you to change the orientation from landscape to portrait to adapt better to the user needs.
- Mount monitors from 19" to 32" with one yoke.
- OASYS monitor yokes offer flexibility through its telescopic abilities and rotation.
- Generous cable raceway allows for multiple cables to be run through the yoke.
- Multiple adaptor options allows connection to OASYS 2001 and OASYS 3001 Spring Arms as well as other arm manufacturers.



Customized OEM solutions available,
contact us for a quote!

Committed to Quality. Committed to You.

SPECIFICATIONS

Monitor sizes	Single - Adjustable up to 32" diagonal, Dual - up to 20" 4:3 monitors
Bolt Pattern	Vesa 75 x 75/100 x 100/200 x 200mm
Brake	Adjustable friction brake
Yoke Rotation	320°, Optional 360° rotation is available with purchase of OASYS multipole slip ring
Monitor Pitch	Customizable up to 320°
Handle	Removable and Sterilizable
Standard Color Options	RAL 9010 / 9002 / 9003 / 7035

